


DISTRIBUIDOR
AUTORIZADO:



netSteril
info@netsteril.com



Design and development of filling machines and innovative packaging solutions



YOUR PRECISION PARTNER
FOR LIQUID FILLING


SPC Engineering,
a subsidiary company of





SEMI-AUTOMATIC FILLING AND STOPPERING MACHINE FOR PREFILLED SYRINGES AND VIALS

PHARMACEUTICALS
R&D
CLINICAL TRIALS
BIOPHARMA
ADVANCED THERAPIES



Maxys[®] semi-automatic machine is used to fill and stopper small-scale vials and syringes productions (clinical batches, R&D etc.).

FEATURES

- Benchtop machine
- One of the most compact machines on the market
- Easily integrated in LAF, RABS or isolators
- H₂O₂ compatible
- Suitable to various Nest formats
- Filling thanks to volumetric piston pump or peristaltic pump
- Filling volume from 0.001 ml to 150 ml
- Output of up to 700 syringes/hr
- CFR 21 Part 11 compliant

DISTRIBUIDOR
AUTORIZADO:



netSteril
info@netsteril.com

104 rue de Lisbonne, Parc d'Activités de SERMENAZ
69 140 RILLIEUX LA PAPE - FRANCE

T: +33 (0)4 78 88 79 56 • F: +33 (0)4 78 88 75 72 • spc@spc.fr • www.spc.fr

