



Manufacturing of high precision tooling for liquid filling

DISTRIBUIDOR
AUTORIZADO:




netSteril
info@netsteril.com



SPC Manufacturing,
a subsidiary company of  SPC
GROUP

YOUR PRECISION PARTNER
FOR LIQUID FILLING



FILLING NEEDLES

PHARMACEUTICALS
R&D
CLINICAL TRIALS
BIOPHARMA
ADVANCED THERAPIES

Our filling, gazing or washing needles are all custom-designed and application dedicated. All models are adaptable to all machines used in pharmaceutical industry.

FEATURES

- 316 L stainless steel design, to meet pharmaceutical requirements
- Internal and external electrolytic polishing (on request), for an optimal surface ($R_a < 0,4$)
- Monobloc design and/or welded design (laser or TIG)
- Decontamination and passivation
- All kinds of tips available (simple-eye, anti-drop, limited tip, spherical tip)
- CGMP compliant
- Traceability documentation supplied
- Possible repairing (up to 50% savings)

DISTRIBUIDOR
AUTORIZADO:



netSteril
info@netsteril.com

104 rue de Lisbonne, Parc d'Activités de SERMENAZ
69 140 RILLIEUX LA PAPE - FRANCE

T : +33 (0)4 78 88 79 56 • F : +33 (0)4 78 88 75 72 • spc@spc.fr • www.spc.fr

