

Design and development of filling machines and innovative packaging solutions



DISTRIBUIDOR
AUTORIZADO



netSteril
info@netsteril.com

SPC Engineering,
a subsidiary company of  SPC GROUP

YOUR PRECISION PARTNER
FOR LIQUID FILLING





HIGH PRECISION NUMERICAL FILLING UNIT

PHARMACEUTICALS
R&D
CLINICAL TRIALS
BIOPHARMA
ADVANCED THERAPIES

MALYS® VOLUMETRIC
PISTON PUMP





With a repeatability rate of +/- 0.1%, our Malys® numerical filling system aims to optimize filling conditions.

FEATURES

- Stand-alone and/or production line upgrade
- Compact size
- Step-by-step or brushless motorization
- Closed loop positioning system
- Aseptic conditions
- Independent adjustment of each pump and of volume by weight
- Precise control of dose ejection, without systematic use of suckback
- Filling volume from 0.001 to 150 mL
- Rate up to 6,000 vials/hr maxi
- Minimal maintenance (for a high usage rate)
- Quick fitting of the pump with the patented system KLYX®
- CFR 21 Part 11 compliant
- IPC option

MALYS® VOLUMETRIC
PERISTALTIC PUMP



DISTRIBUIDOR
AUTORIZADO



netSteril
info@netsteril.com

104 rue de Lisbonne, Parc d'Activités de SERMENAZ
69 140 RILLIEUX LA PAPE - FRANCE

T : +33 (0)4 78 88 79 56 • F : +33 (0)4 78 88 75 72 • spc@spc.fr • www.spc.fr

